



Testing of Peroxide Forming Chemicals Using MQuant Test Strips

1. Visually inspect bottle for:
 - a. Crystallization (especially around cap or in the liquid)
 - b. Visible discoloration
 - c. Liquid stratification
2. Use glass Pasteur pipet to get small amount of liquid
3. Put drop on square test field on test strip
4. Shake off excess liquid
5. Add one drop of DI water when testing an organic solvent
 - a. (Do not need to add any water for inorganic solutions)
6. Write results on label
 - a. Include date, value, and your initials
7. Allow test strip to evaporate, then throw away in the trash

QUANTITATIVE TESTING

AFTER VISUAL INSPECTION, COMPLETE QUANTITATIVE TESTING. IF PEROXIDES ARE DETECTED, FOLLOW THESE RECOMMENDATIONS FOR CONTINUED USE OR DISPOSAL:

1 PPM	Okay to keep; Do not concentrate to near-dryness.
3 PPM	Okay to keep; Do not concentrate to near-dryness.
10 PPM	Okay to keep, but must now test every month; Do not concentrate to near-dryness.
>10 PPM	Get 0-25ppm test strips from Stockroom to get accurate readings; Do not concentrate to near-dryness.
≥20 PPM	19ppm is the maximum level that the hazardous waste hauler will transport. Combine, bulk, or inhibit/stabilize until ppm is 19ppm or less. Stockroom has test strips that read more accurately between 0-25 ppm. Do not concentrate to near-dryness.
>20 PPM <100 PPM	Not recommended for distilling or concentrating.
≥100 PPM	DO NOT USE. Contact EHS for assistance. Additional testing may be required prior to disposal.